

PE anti-human CD274 (B7-H1, PD-L1)

Mouse IgG2b, kappa
29E.2A3

BioLegend
From Revvity



25 tests

Artikel-Nr.: B329705 | **BioLegend** | Herstellernummer: 329705

71,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC

Clone: 29E.2A3

Format: PE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 tests

PE anti-human CD274 (B7-H1, PD-L1)

Clone: 29E.2A3

Isotype: Mouse IgG2b, kappa

Other_Names: Programmed cell death ligand 1 (PD-L1), B7 homolog 1 (B7-H1)

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=4375>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>