

96er DeepWell Platte riplate medio, 1.2 ml, RW

U-Boden, Polypropylen
SBS Format

Packung à 32 Stück



Artikel-Nr.: 742917 | **Ritter** | Herstellernummer: 43001-0101

154,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Boden: U-Boden

Eigenschaft: Medio

Kategorie: 96er DeepWell Platte

Material: PP

Verpackung: Packung à 32 Stück

Volumen: 0.5 ml

Well Form: rund

Riplate® medio

Bis zu 50% weniger Stapelhöhe

- Durch sich tangential berührende Wells benötigt Riplate® medio bis zu 50% geringere Stapelhöhe.
- Die erhöhten Ränder der einzelnen Wells lassen die Riplate® medio luftdicht verschweißen oder verkleben.
- Die alphanumerische Kennzeichnung ermöglicht die Registrierung und Identifizierung der einzelnen Proben.
- Neukonzipierte Wells ermöglichen eine gute Probendurchmischung.

Riplate® medio 1.0 ml

- Höhe: nur 27.0 mm
- Gesamtvolumen: 1.2 ml / Well
- Arbeitsvolumen: 50 – 1000 µl

Riplate® medio

- *Due to the tangency of Riplate® medio, less than 50 % stacking size is needed.*
- *Higher well rims allow an airtight closure or sealing of Riplate® medio.*
- *The alphanumeric marking allows the registration and identification of single samples.*
- *New developed conical wells provide a good bioturbation.*

Riplate® medio 1.0 ml

- *Height: only 27.0 mm*
- *Total volume: 1.2 ml / well*
- *Working volume: 50 - 1000 µl / well*