

PE anti-human CD167a (DDR1)

Mouse IgM, kappa
51D6

BioLegend
From Revvity



100 tests

Artikel-Nr.: B334006 | **BioLegend** | Herstellernummer: 334006

175,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC

Clone: 51D6

Format: PE

Isotype: kappa, Mouse IgM

Species of Origin: Mouse

Verpackung: 100 Tests

PE anti-human CD167a (DDR1)

Clone: 51D6

Isotype: Mouse IgM, kappa

Other_Names: Discoidin domain receptor, DDR1, neuroepithelial tyrosine kinase, epithelial-specific receptor kinase, tyrosine kinase receptor E (trkE), cell adhesion kinase (CAK), PTK3, RTK6

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=4903>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>