

## APC anti-human C3AR

Mouse IgG2b, kappa  
hC3aRZ8

**BioLegend**  
From Revvity



25 tests

Artikel-Nr.: B345805 | BioLegend | Herstellernummer: 345805

**124,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FC  
**Clone:** hC3aRZ8  
**Format:** APC  
**Isotype:** kappa, Mouse IgG2b  
**Species of Origin:** Mouse  
**Verpackung:** 25 tests

APC anti-human C3AR

Clone: hC3aRZ8  
Isotype: Mouse IgG2b, kappa  
Other\_Names: C3a anaphylatoxin chemotactic receptor, C3a-R, AZ3B, HNFAG09  
Size: 25 tests

Reactivity: Human  
Format: APC  
App\_Abbrev: FC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=7274>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>