

PE anti-human CD85d (ILT4)

Rat IgG2a, kappa
42D1

BioLegend
From Revvity



100 tests

Artikel-Nr.: B338706 | BioLegend | Herstellernummer: 338706

254,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: 42D1

Format: PE

Isotype: kappa, Rat IgG2a

Species of Origin: Rat

Verpackung: 100 Tests

PE anti-human CD85d (ILT4)

Clone: 42D1

Isotype: Rat IgG2a, kappa

Other_Names: CD85, ILT4, MIR10, Leukocyte Immunoglobulin-Like Receptor subfamily B member 2 (LILRB2), Leukocyte Immunoglobulin-like Receptor 2 (LIR-2)

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=5576>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>