

PerCP/Cyanine5.5 anti-human CD271 (NGFR)

Mouse IgG1, kappa
ME20.4

BioLegend
From Revvity



100 tests

Artikel-Nr.: B345112 | [BioLegend](#) | Herstellernummer: 345112

306,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: ME20.4
Format: PERCP_CYANINE5.5
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 Tests

PerCP/Cyanine5.5 anti-human CD271 (NGFR)

Clone: ME20.4
Isotype: Mouse IgG1, kappa
Other_Names: p75NTR, TNFRSF16, p75(NTR), Gp80-LNGFR, NGFR, Nerve Growth Factor Receptor
Size: 100 tests

Reactivity: Human
Format: PERCP_CYANINE5.5
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=7834>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>