

PE anti-human Ksp37 (FGFBP2)

Mouse IgG1, kappa
TDA3

BioLegend
From Revvity



25 tests

Artikel-Nr.: B346603 | **BioLegend** | Herstellernummer: 346603

87,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICFC

Clone: TDA3

Format: PE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 25 tests

PE anti-human Ksp37 (FGFBP2)

Clone: TDA3

Isotype: Mouse IgG1, kappa

Other_Names: Killer-specific secretory protein of 37 kD (Ksp37), FGF binding protein 2 (FGFBP2), Fibroblast growth factor binding protein 2

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6328>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>