

Purified anti-human CD254 (TRANCE, RANKL)

Mouse IgG2b, kappa
MIH24

BioLegend
From Revvity



25 µg

Artikel-Nr.: B347501 | **BioLegend** | Herstellernummer: 347501

100,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA, FC, IHC-P

Clone: MIH24

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-human CD254 (TRANCE, RANKL)

Clone: MIH24

Isotype: Mouse IgG2b, kappa

Other_Names: TRANCE (TNF-related activation-induced cytokine receptor), CD254, RANK-L (receptor activator of NF-kappaB ligand), OPGL (Osteoprotegerin ligand), SOFA, TNFSF-11

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ELISA, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6396>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>