

LEGEND MAX Human anti-mouse Ig (HAMA) ELISA Kit

BioLegend
From Revvity



1 Pre-coated Plate

Artikel-Nr.: B438307 | BioLegend | Herstellernummer: 438307

246,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ELISA
Format: LEGEND MAX
Species of Origin: N/A
Verpackung: 1 Pre-coated Plate

LEGEND MAX Human anti-mouse Ig (HAMA) ELISA Kit

Clone:
Isotype:
Other_Names: Human anti-mouse antibodies (HAMA) Pre-coated ELISA Kit
Size: 1 Pre-coated Plate

Reactivity: Human
Format: LEGEND MAX
App_Abbrev: ELISA
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=7266>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>