

Purified anti-human CD134 (OX40)

Mouse IgG1, kappa
Ber-ACT35 (ACT35)

BioLegend
From Revvity



100 µg

Artikel-Nr.: B350002 | **BioLegend** | Herstellernummer: 350002

71,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ELISA, FC, IHC-F, IHC-P, WB

Clone: Ber-ACT35 (ACT35)

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human CD134 (OX40)

Clone: Ber-ACT35 (ACT35)

Isotype: Mouse IgG1, kappa

Other_Names: ACT35 antigen, tumor necrosis factor receptor superfamily member 4 (TNFRSF4)

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, IHC-P, ELISA, IHC-F, WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6802>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>