

FITC anti-human MRP-14 (S100A9)

Mouse IgG1, kappa
MRP 1H9

BioLegend
From Revvity



100 tests

Artikel-Nr.: B350704 | **BioLegend** | Herstellernummer: 350704

152,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICFC

Clone: MRP 1H9

Format: FITC

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 Tests

FITC anti-human MRP-14 (S100A9)

Clone: MRP 1H9

Isotype: Mouse IgG1, kappa

Other_Names: S100A9, Migration inhibitory factor-related proteins-14, Calgranulin-B, Leukocyte L1 complex heavy chain, Calprotectin L1H subunit

Size: 100 tests

Reactivity: Human

Format: FITC

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6977>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>