

APC anti-human CD201 (EPCR)

Rat IgG1, kappa

RCR-401

BioLegend
From Revvity



25 tests

Artikel-Nr.: B351905 | [BioLegend](#) | Herstellernummer: 351905

91,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: RCR-401

Format: APC

Isotype: kappa, Rat IgG1

Species of Origin: Rat

Verpackung: 25 tests

APC anti-human CD201 (EPCR)

Clone: RCR-401

Isotype: Rat IgG1, kappa

Other_Names: Endothelial protein C receptor, protein C receptor, PROCR, activated protein C receptor, APC receptor

Size: 25 tests

Reactivity: Human

Format: APC

App_Abbrev: FC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=7693>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>