

Purified anti-mouse/human Mac-2 (Galectin-3)

Mouse IgG1, kappa
Gal397

BioLegend
From Revvity



500 µg

Artikel-Nr.: B126702 | **BioLegend** | Herstellernummer: 126702

155,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ICC, WB
Clone: Gal397
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-mouse/human Mac-2 (Galectin-3)

Clone: Gal397
Isotype: Mouse IgG1, kappa
Other_Names: Galectin-3, Mac-2, Gal-3, RL-29, galactose-specific lectin 3, CBP-35, L-34, epsilonBP
Size: 500 µg

Reactivity: Mouse, Human
Format: PURE
App_Abbrev: WB, ICC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=4900>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>