

Purified anti-human IFN-beta

Mouse IgG1, kappa
IFNb/A1

BioLegend
From Revvity



200 µg

Artikel-Nr.: B514002 | BioLegend | Herstellernummer: 514002

235,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA, Neut, WB

Clone: IFNb/A1

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 200 µg

Purified anti-human IFN-beta

Clone: IFNb/A1

Isotype: Mouse IgG1, kappa

Other_Names: Interferon-beta, Fibroblast interferon (Fi-IFN, F-IFN), Type-1 interferon, acid-stable interferon, R1-GI factor, IFN-b, IFN-beta

Size: 200 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, ELISA, Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=5133>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>