

Purified anti-mouse CD279 (PD-1)

Rat IgG2a, kappa
29F.1A12

BioLegend
From Revvity



500 µg

Artikel-Nr.: B135202 | BioLegend | Herstellernummer: 135202

109,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Block, ELISA, FC, IHC-F

Clone: 29F.1A12

Format: PURE

Isotype: kappa, Rat IgG2a

Species of Origin: Rat

Verpackung: 500 µg

Purified anti-mouse CD279 (PD-1)

Clone: 29F.1A12

Isotype: Rat IgG2a, kappa

Other_Names: PD-1, Programmed Death-1, PDCD1

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, IHC-F, Block, ELISA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6167>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>