

PE anti-human CD172a/b (SIRPalpha/beta)

Mouse IgG1, kappa
SE5A5

BioLegend
From Revvity



100 tests

Artikel-Nr.: B323806 | **BioLegend** | Herstellernummer: 323806

139,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC

Clone: SE5A5

Format: PE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 Tests

PE anti-human CD172a/b (SIRPalpha/beta)

Clone: SE5A5

Isotype: Mouse IgG1, kappa

Other_Names: SIRPa, SIRPb, SIRPalpha/beta, BIT, SHPS1, MFR, P84, PTPNS1, CD172 antigen-like family member A, CD172 antigen-like family member B

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=4030>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>