

Ultra-LEAF Purified anti-human HLA-A,B,C

Mouse IgG2a, kappa
W6/32

BioLegend
From Revvity



1 mg

Artikel-Nr.: B311428 | **BioLegend** | Herstellernummer: 311428

90,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Activ, Block, FC, IHC-F, IP, WB

Clone: W6/32

Format: ULEAF

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 1 mg

Ultra-LEAF Purified anti-human HLA-A,B,C

Clone: W6/32

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Class I, MHC class I

Size: 1 mg

Reactivity: Human, Cynomolgus, Rhesus

Format: ULEAF

App_Abbrev: FC, Activ, Block, IHC-F, IP, WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8095>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>