

Purified anti-human EphA2

Mouse IgG2b, kappa
SHM16

BioLegend
From Revvity



25 µg

Artikel-Nr.: B356801 | [BioLegend](#) | Herstellernummer: 356801

67,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC

Clone: SHM16

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-human EphA2

Clone: SHM16

Isotype: Mouse IgG2b, kappa

Other_Names: Ephrin type-A receptor 2 (EPHA2), Epithelial cell kinase, Tyrosine-protein kinase receptor ECK, Myk2, Sek2, EPH receptor A2, eEpithelial cell receptor protein tyrosine kinase, Protein tyrosine kinase, Soluble EPHA2 variant 1

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8365>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>