

## PE anti-human EphA2

Mouse IgG2b, kappa  
SHM16

**BioLegend**  
From Revvity



25 tests

Artikel-Nr.: B356803 | **BioLegend** | Herstellernummer: 356803

**87,75 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Aktionsangebote:** Aktion

**App\_Abbrev:** FC

**Clone:** SHM16

**Format:** PE

**Isotype:** kappa, Mouse IgG2b

**Species of Origin:** Mouse

**Verpackung:** 25 tests

PE anti-human EphA2

Clone: SHM16

Isotype: Mouse IgG2b, kappa

Other\_Names: Ephrin type-A receptor 2 (EPHA2), Epithelial cell kinase, Tyrosine-protein kinase receptor ECK, Myk2, Sek2, EPH receptor A2, eEpithelial cell receptor protein tyrosine kinase, Protein tyrosine kinase, Soluble EPHA2 variant 1

Size: 25 tests

Reactivity: Human

Format: PE

App\_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8366>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>