

## Purified anti-Rad23B

Mouse IgG1, lambda  
11A2SC

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B654201 | **BioLegend** | Herstellernummer: 654201

**118,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** KO/KD-WB, WB

**Clone:** 11A2SC

**Format:** PURE

**Isotype:** lambda, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 25 µg

Purified anti-Rad23B

Clone: 11A2SC

Isotype: Mouse IgG1, lambda

Other\_Names: RAD23 homolog B (S. cerevisiae), RAD23, yeast homolog of B, UV excision repair protein RAD23 homolog B, HHR23B

Size: 25 µg

Reactivity: Mouse, Human

Format: PURE

App\_Abbrev: WB, KO/KD-WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8355>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>