

Purified anti-TCF1 (TCF7)

Mouse IgG1, kappa
7F11A10

BioLegend
From Revvity



100 µg

Artikel-Nr.: B655202 | BioLegend | Herstellernummer: 655202

231,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, IHC-P, WB

Clone: 7F11A10

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-TCF1 (TCF7)

Clone: 7F11A10

Isotype: Mouse IgG1, kappa

Other_Names: T-cell-specific transcription factor 1, transcription factor 7, TCF7

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, ICC, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8541>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>