

## Purified anti-Tubulin Beta 3 (TUBB3)

Mouse IgG2a, kappa  
AA10

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B657401 | **BioLegend** | Herstellernummer: 657401

**60,75 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** ICC, WB

**Clone:** AA10

**Format:** PURE

**Isotype:** kappa, Mouse IgG2a

**Species of Origin:** Mouse

**Verpackung:** 25 µg

Purified anti-Tubulin Beta 3 (TUBB3)

Clone: AA10

Isotype: Mouse IgG2a, kappa

Other\_Names: beta-3 Tubulin, Neuronal class III beta-tubulin, beta3-tub

Size: 25 µg

Reactivity: Mouse, Rat, Human

Format: PURE

App\_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8921>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>