

PerCP/Cyanine5.5 anti-human RANTES (CCL5)

Mouse IgG2b, kappa
VL1

BioLegend
From Revvity



100 tests

Artikel-Nr.: B515508 | **BioLegend** | Herstellernummer: 515508

306,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC
Clone: VL1
Format: PERCP_CYANINE5.5
Isotype: kappa, Mouse IgG2b
Species of Origin: Mouse
Verpackung: 100 Tests

PerCP/Cyanine5.5 anti-human RANTES (CCL5)

Clone: VL1
Isotype: Mouse IgG2b, kappa
Other_Names: CCL5, C-C motif chemokine 5, Small-inducible cytokine A5, T-cell-specific protein RANTES, SIS-delta, T cell-specific protein P228, TCP228, Eosinophil-chemotactic cytokine, EoCP
Size: 100 tests

Reactivity: Human
Format: PERCP_CYANINE5.5
App_Abbrev: ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=8981>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>