

Purified anti-human CCL17 (TARC)

Goat Polyclonal Ig

Poly5230

BioLegend
From Revvity



100 µg

Artikel-Nr.: B523002 | BioLegend | Herstellernummer: 523002

280,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Capture, WB

Clone: Poly5230

Format: PURE

Isotype: Goat Polyclonal Ig

Species of Origin: Goat

Verpackung: 100 µg

Purified anti-human CCL17 (TARC)

Clone: Poly5230

Isotype: Goat Polyclonal Ig

Other_Names: Thymus and activation-regulated chemokine (TARC), Small-inducible cytokine 17 (SCYA17), ABCD-2

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture, WB

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9156>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>