

Biotin anti-human CCL17 (TARC)

Goat Polyclonal Ig

Poly5230

BioLegend
From Revvity



50 µg

Artikel-Nr.: B523003 | BioLegend | Herstellernummer: 523003

300,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Detection

Clone: Poly5230

Format: BIOTIN

Isotype: Goat Polyclonal Ig

Species of Origin: Goat

Verpackung: 50 µg

Biotin anti-human CCL17 (TARC)

Clone: Poly5230

Isotype: Goat Polyclonal Ig

Other_Names: Thymus and activation-regulated chemokine (TARC), Small-inducible cytokine 17 (SCYA17), ABCD-2

Size: 50 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: ELISA Detection

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9169>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>