

Alexa Fluor® 488 anti-Tubulin Beta 3 (TUBB3)

Mouse IgG2a, kappa
AA10

BioLegend
From Revvity



25 µg

Artikel-Nr.: B657403 | **BioLegend** | Herstellernummer: 657403

105,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: 3D IHC, ICFC, IHC-F

Clone: AA10

Format: A488

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 25 µg

Alexa Fluor® 488 anti-Tubulin Beta 3 (TUBB3)

Clone: AA10

Isotype: Mouse IgG2a, kappa

Other_Names: beta-3 Tubulin, Neuronal class III beta-tubulin, beta3-tub

Size: 25 µg

Reactivity: Mouse, Rat, Human

Format: A488

App_Abbrev: ICFC, IHC-F, 3D IHC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9001>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>