

Purified anti-mouse CD25 (Maxpar® Ready)

Rat IgG2b, kappa
3C7

BioLegend
From Revvity



100 µg

Artikel-Nr.: B101913 | **BioLegend** | Herstellernummer: 101913

95,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: CyTOF®, FC

Clone: 3C7

Format: PUREMAX

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 100 µg

Purified anti-mouse CD25 (Maxpar® Ready)

Clone: 3C7

Isotype: Rat IgG2b, kappa

Other_Names: IL-2Ralpha, Ly-43, p55, Tac

Size: 100 µg

Reactivity: Mouse

Format: PUREMAX

App_Abbrev: FC, CyTOF®

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9160>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>