

PE anti-human RANTES (CCL5)

Mouse IgG2b, kappa
VL1

BioLegend
From Revvity



25 tests

Artikel-Nr.: B515503 | BioLegend | Herstellernummer: 515503

104,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC

Clone: VL1

Format: PE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 tests

PE anti-human RANTES (CCL5)

Clone: VL1

Isotype: Mouse IgG2b, kappa

Other_Names: CCL5, C-C motif chemokine 5, Small-inducible cytokine A5, T-cell-specific protein RANTES, SIS-delta, T cell-specific protein P228, TCP228, Eosinophil-chemotactic cytokine, EoCP

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=6112>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>