

APC anti-human/mouse/rat MR1

Mouse IgG2a, kappa
26.5

BioLegend
From Revvity



25 tests

Artikel-Nr.: B361107 | [BioLegend](#) | Herstellernummer: 361107

124,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: 26.5

Format: APC

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 25 tests

APC anti-human/mouse/rat MR1

Clone: 26.5

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Complex Class I-related, HLALS

Size: 25 tests

Reactivity: Human, Mouse, Rat

Format: APC

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9373>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>