

Purified anti-human CD130 (gp130)

Mouse IgG2a, kappa
2E1B02

BioLegend
From Revvity



100 µg

Artikel-Nr.: B362002 | **BioLegend** | Herstellernummer: 362002

100,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: 2E1B02

Format: PURE

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human CD130 (gp130)

Clone: 2E1B02

Isotype: Mouse IgG2a, kappa

Other_Names: GP130, GP-130, IL6ST, IL-6Rbeta, Oncostatin M Receptor, Oncostatin-M receptor subunit alpha

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9377>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>