

LEGEND MAX Human Adiponectin (Acrp30) ELISA Kit

BioLegend
From Revvity



1 Pre-coated Plate

Artikel-Nr.: B442307 | [BioLegend](#) | Herstellernummer: 442307

398,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA
Format: LEGEND MAX
Species of Origin: N/A
Verpackung: 1 Pre-coated Plate

LEGEND MAX Human Adiponectin (Acrp30) ELISA Kit

Clone:

Isotype:

Other_Names: Adiponectin, Adipocyte complement-related protein of 30 kD (Acrp30), Adipose most abundant gene transcript 1 (apM1), Gelatin-binding protein of 28 kD (GBP28), ADIPOQ

Size: 1 Pre-coated Plate

Reactivity: Human
Format: LEGEND MAX
App_Abbrev: ELISA
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9410>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>