

## Purified anti-NLRC4 (IPAF)

Mouse IgG2a, kappa  
6H9B13

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B659702 | **BioLegend** | Herstellernummer: 659702

**159,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**Aktionsangebote:** Aktion

**App\_Abbrev:** ICC, WB

**Clone:** 6H9B13

**Format:** PURE

**Isotype:** kappa, Mouse IgG2a

**Species of Origin:** Mouse

**Verpackung:** 100 µg

Purified anti-NLRC4 (IPAF)

Clone: 6H9B13

Isotype: Mouse IgG2a, kappa

Other\_Names: NLR family CARD domain-containing protein 4, CARD, LRR, and NACHT-containing protein (CLAN), Caspase recruitment domain-containing protein 12 (CARD12), Ice protease-activating factor (IPAF)

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9599>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>