

Purified anti-CCL11 (Eotaxin)

Rat IgG2a, kappa
L402H11

BioLegend
From Revvity



100 µg

Artikel-Nr.: B660402 | **BioLegend** | Herstellernummer: 660402

203,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: WB
Clone: L402H11
Format: PURE
Isotype: kappa, Rat IgG2a
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-CCL11 (Eotaxin)

Clone: L402H11
Isotype: Rat IgG2a, kappa
Other_Names: Chemokine (C-C Motif) ligand 11, Eotaxin-1, Eosinophil chemotactic protein, Small-inducible cytokine A11
Size: 100 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: WB
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=9723>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>