


PCR TW 12-Tube Strip, 0.2 ml, frosted

Extra Robust, High Profile

DNA-, DNase-, RNA-, RNase-, Pyrogenfrei



Beutel à 80 Strips

Artikel-Nr.: 712183 | Biozym

91,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Farbe: frosted

Kategorie: 12er-Streifen

Produktlinie: Premium Line

Profil: High Profile

Verpackung: Beutel à 80 Strips

Volumen: 0.2 ml