

## Recombinant Mouse LIF (carrier-free)

**BioLegend**  
From Revvity



10 µg

Artikel-Nr.: B554002 | [BioLegend](#) | Herstellernummer: 554002

**113,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

App\_Abbrev: BA

Format: RECOM

Verpackung: 10 µg

Recombinant Mouse LIF (carrier-free)

Clone:

Isotype:

Other\_Names: MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor), LIF (leukemia inhibitory factor)

Size: 10 µg

Reactivity: Mouse

Format: RECOM

App\_Abbrev: BA

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10032>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>