

Biotin anti-human CCL5 (RANTES)

Goat Polyclonal Ig
Poly5197

BioLegend
From Revvity



50 µg

Artikel-Nr.: B519703 | BioLegend | Herstellernummer: 519703

300,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Detection

Clone: Poly5197

Format: BIOTIN

Isotype: Goat Polyclonal Ig

Species of Origin: Goat

Verpackung: 50 µg

Biotin anti-human CCL5 (RANTES)

Clone: Poly5197

Isotype: Goat Polyclonal Ig

Other_Names: SIS-delta, T-cell specific protein p288, beta-chemokine RANTES, regulated upon activation, normally T-expressed, and presumably secreted, small inducible cytokine A5

Size: 50 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: ELISA Detection

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10633>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>