

## APC anti-human CD19

Mouse IgG1, kappa  
SJ25C1

**BioLegend**  
From Revvity



25 tests

Artikel-Nr.: B363005 | **BioLegend** | Herstellernummer: 363005

**80,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FC  
**Clone:** SJ25C1  
**Format:** APC  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 25 tests

APC anti-human CD19

Clone: SJ25C1  
Isotype: Mouse IgG1, kappa  
Other\_Names: CD19 Molecule, CD19 Antigen, Differentiation Antigen CD19, B-Lymphocytes Surface Antigen B4, T-cell Surface Antigen Leu-12, CVID3, B4, B-lymphocyte Antigen CD19  
Size: 25 tests

Reactivity: Human  
Format: APC  
App\_Abbrev: FC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10264>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>