

Purified anti-Nanog

Mouse IgG1, kappa
23D2-3C6

BioLegend
From Revvity



25 µg

Artikel-Nr.: B674001 | **BioLegend** | Herstellernummer: 674001

88,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ICC, ICFC, WB
Clone: 23D2-3C6
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Nanog

Clone: 23D2-3C6
Isotype: Mouse IgG1, kappa
Other_Names: Homeobox Transcription Factor Nanog, Nanog Homeobox, Homeobox Protein NANOG, Homeobox Transcription Factor Nanog-Delta 48
Size: 25 µg

Reactivity: Human
Format: PURE
App_Abbrev: WB, ICFC, ICC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10616>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>