

## Purified anti-Nanog

Mouse IgG1, kappa  
23D2-3C6

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B674002 | **BioLegend** | Herstellernummer: 674002

**186,75 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**Aktionsangebote:** Aktion  
**App\_Abbrev:** ICC, ICFC, WB  
**Clone:** 23D2-3C6  
**Format:** PURE  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-Nanog

Clone: 23D2-3C6  
Isotype: Mouse IgG1, kappa  
Other\_Names: Homeobox Transcription Factor Nanog, Nanog Homeobox, Homeobox Protein NANOG, Homeobox Transcription Factor Nanog-Delta 48  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: WB, ICFC, ICC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10616>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>