

PerCP/Cyanine5.5 anti-human TSLPR (TSLP-R)

Mouse IgG1, kappa
1B4

BioLegend
From Revvity



100 tests

Artikel-Nr.: B322812 | **BioLegend** | Herstellernummer: 322812

306,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: 1B4
Format: PERCP_CYANINE5.5
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 Tests

PerCP/Cyanine5.5 anti-human TSLPR (TSLP-R)

Clone: 1B4
Isotype: Mouse IgG1, kappa
Other_Names: Thymic stromal lymphopoietin protein receptor, cytokine receptor-like 2, CRL2, CRLF2, IL-XR, TSLP-R
Size: 100 tests

Reactivity: Human
Format: PERCP_CYANINE5.5
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10409>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>