

Purified anti-LIF

Rat IgG2a, kappa
M1506B09

BioLegend
From Revvity



100 µg

Artikel-Nr.: B674702 | **BioLegend** | Herstellernummer: 674702

138,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICC, WB

Clone: M1506B09

Format: PURE

Isotype: kappa, Rat IgG2a

Species of Origin: Rat

Verpackung: 100 µg

Purified anti-LIF

Clone: M1506B09

Isotype: Rat IgG2a, kappa

Other_Names: Leukemia inhibitory factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 100 µg

Reactivity: Mouse, Human

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=10856>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>