

## Purified anti-Aggrecan

Mouse IgG1, kappa  
M4004A03

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B674902 | BioLegend | Herstellernummer: 674902

**212,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** WB  
**Clone:** M4004A03  
**Format:** PURE  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-Aggrecan

Clone: M4004A03  
Isotype: Mouse IgG1, kappa  
Other\_Names: Aggrecan 1, AGC1, Aggrecan Proteoglycan, Chondroitin Sulfate Proteoglycan 1, Cartilage-Specific Proteoglycan Core Protein, CSPG1, MSK16, CSPGCP, SEDK, and CSPCP  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: WB  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11665>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>