

Purified anti-CDKN2A (p16)

Mouse IgG2b, kappa
15C10C30

BioLegend
From Revvity



100 µg

Artikel-Nr.: B675602 | **BioLegend** | Herstellernummer: 675602

249,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, KO/KD-WB, WB

Clone: 15C10C30

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-CDKN2A (p16)

Clone: 15C10C30

Isotype: Mouse IgG2b, kappa

Other_Names: Cyclin-dependent kinase inhibitor 2A, Cyclin-dependent kinase 4 inhibitor A (CDK4I), Multiple tumor suppressor 1 (MTS-1), p16-INK4a, p14 ARF

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, ICC, KO/KD-WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11761>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>