

## Anti-CD230 (Prion)

Mouse IgG2a, kappa  
3F4

**BioLegend**  
From Revvity



200 µL

Artikel-Nr.: B800307 | [BioLegend](#) | Herstellernummer: 800307

**844,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** FC

**Clone:** 3F4

**Format:** ASCITES

**Isotype:** kappa, Mouse IgG2a

**Species of Origin:** Mouse

**Verpackung:** 200 µl

Anti-CD230 (Prion)

Clone: 3F4

Isotype: Mouse IgG2a, kappa

Other\_Names: Major prion protein, p27-30, CD230 antigen, prion protein PrP, prion-related protein

Size: 200 µL

Reactivity: Human

Format: ASCITES

App\_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11027>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>