

Purified (azide-free) anti-beta-Amyloid, 1-16

Mouse IgG1, kappa
6E10

BioLegend
From Revvity



200 µL

Artikel-Nr.: B803001 | **BioLegend** | Herstellernummer: 803001

444,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA, EM, ICC, IHC-F, IHC-P, WB

Clone: 6E10

Format: PUREAF

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 200 µl

Purified (azide-free) anti-beta-Amyloid, 1-16

Clone: 6E10

Isotype: Mouse IgG1, kappa

Other_Names: AAA, ABETA, ABPP, AD1, APPI, CTFgamma, CVAP, PN-II, PN2, Amyloid beta A4 protein, preA4, protease, peptidase nexin-II, beta-amyloid peptide, alzheimer disease amyloid protein, cerebral vascular amyloid peptide, APP, Amyloid Precursor Protein

Size: 200 µL

Reactivity: Human

Format: PUREAF

App_Abbrev: WB, Direct ELISA, IHC-P, IHC-F, EM, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11228>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>