

Purified (azide-free) anti-Tau, x-421

Mouse IgG1
C3

BioLegend
From Revvity



200 µL

Artikel-Nr.: B806301 | [BioLegend](#) | Herstellernummer: 806301

244,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Direct ELISA, IHC-P, WB

Clone: C3

Format: PUREAF

Isotype: Mouse IgG1

Species of Origin: Mouse

Verpackung: 200 µl

Purified (azide-free) anti-Tau, x-421

Clone: C3

Isotype: Mouse IgG1

Other_Names: Microtubule-associated protein tau, PHF-tau, paired helical filament-tau, neurofibrillary tangle, microtubule-associated protein tau, isoform 4, G protein beta1/gamma2 subunit-interacting factor 1, DDPAC, FTDP-17, MAPTL, MSTD, MTBT1, MTBT2, PPND
Size: 200 µL

Reactivity: Human

Format: PUREAF

App_Abbrev: IHC-P, WB, Direct ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11570>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>