

Purified (azide-free) anti-AICD

Rabbit Polyclonal IgG

Poly8119

BioLegend
From Revvity



50 µg

Artikel-Nr.: B811901 | **BioLegend** | Herstellernummer: 811901

317,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ELISA

Clone: Poly8119

Format: PUREAF

Isotype: Rabbit Polyclonal IgG

Species of Origin: Rabbit

Verpackung: 50 µg

Purified (azide-free) anti-AICD

Clone: Poly8119

Isotype: Rabbit Polyclonal IgG

Other_Names: Amyloid beta A4 protein, preA4, protease nexin-II, peptidase nexin-II, beta-amyloid peptide, alzheimer disease amyloid protein, cerebral vascular amyloid peptide, APP, Amyloid Precursor Protein

Size: 50 µg

Reactivity: Human

Format: PUREAF

App_Abbrev: ELISA

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11253>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>