

Purified (azide-free) anti-PARIS (ZNF746)

Mouse IgG1, kappa
N196/16

BioLegend
From Revvity



100 µL

Artikel-Nr.: B820601 | **BioLegend** | Herstellernummer: 820601

99,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: IHC-P, WB

Clone: N196/16

Format: PUREAF

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µl

Purified (azide-free) anti-PARIS (ZNF746)

Clone: N196/16

Isotype: Mouse IgG1, kappa

Other_Names: Zinc finger protein 746, parkin-interacting sustrate, parkin-interacting substrate, ZNF746

Size: 100 µL

Reactivity: Human, Mouse, Rat

Format: PUREAF

App_Abbrev: IHC-P, WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11507>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>