

Purified anti-V5 Epitope Tag

Rabbit Polyclonal IgG

Poly29038

BioLegend
From Revvity



100 µL

Artikel-Nr.: B903801 | [BioLegend](#) | Herstellernummer: 903801

265,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA, ICC, IHC, Purification, WB

Clone: Poly29038

Format: PURE

Isotype: Rabbit Polyclonal IgG

Species of Origin: Rabbit

Verpackung: 100 µl

Purified anti-V5 Epitope Tag

Clone: Poly29038

Isotype: Rabbit Polyclonal IgG

Other_Names: GKPIPPLLGLDST epitope tag, GKPIPPLLGLDST tag, Protein Rev, Regulator of expression of viral proteins, Rev, V5 epitope tag

Size: 100 µL

Reactivity:

Format: PURE

App_Abbrev: WB, ICC, IHC, ELISA, Purification

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11593>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>