

Anti-CD31 (PECAM-1)

Mouse IgG1, kappa
P2B1

BioLegend
From Revvity



100 µL

Artikel-Nr.: B910003 | BioLegend | Herstellernummer: 910003

203,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, IHC-F, IP, WB

Clone: P2B1

Format: ASCITES

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µl

Anti-CD31 (PECAM-1)

Clone: P2B1

Isotype: Mouse IgG1, kappa

Other_Names: Platelet endothelial cell adhesion molecule, PECA1, GPIIA', endoCAM, CD31 antigen, CD31/EndoCAM, adhesion molecule, PECAM-1, CD31/EndoCAM

Size: 100 µL

Reactivity: Human

Format: ASCITES

App_Abbrev: IHC-F, ICC, WB, IP

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11135>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>